

Fast DNA-spin™ Plasmid DNA Purification Kit

For the efficient isolation plasmid DNA



PRODUCT FEATURES

- Fast DNA-spin™ Plasmid DNA Purification Kit is designed for rapid purification of plasmids from bacterial cells.
- The fastest, simplest procedure for purifying the highest quality plasmid DNA.
- The high quality plasmid DNA can be used directly for the downstream application.
- Minimal nicks of plasmid DNA guarantees good results in plasmid DNA sequencing.

INTRODUCTION

The method that used in Fast DNA-spin™ Plasmid DNA Purification Kit is optimized from traditional alkaline lysis technology by which high quality plasmid DNA could be purified within 10 minutes. The new lysis buffer allows the adsorption of DNA onto silica membrane in the presence of high salt. The material that is used to make the silica membrane is unique, highly-efficient and highly-specified. This protocol is designed for purification of DNA from 1-4 ml overnight cultures of *E. coli*. Plasmid DNA prepared by Fast DNA-spin™ Plasmid DNA Purification Kit is suitable for a variety of routine applications including restriction enzyme digestion, sequencing, library screening, ligation and transformation, *in vitro* translation and transfection of robust cells.

KIT CONTENTS

Components	50 prep
Resuspension buffer ¹ (Bottle 1)	15 ml
Lysis buffer ² (Bottle 2)	15 ml
Neutralization buffer (Bottle 3)	20 ml
Washing Buffer ³ (Bottle 4)	15 ml
Elution Buffer (Bottle 5)	15 ml
Lysis Viewer	75 µl
RNase A ⁴ (10 mg/ml)	150 µl
Spin Column	50 ea
Collection Tube	50 ea
Storage Conditions : Room Temperature	

¹ Briefly spin the dissolved RNase A solution and add the RNase A solution to Resuspension Buffer. Before use, store Resuspension Buffer at 4 °C after adding RNase A solution.

² Check Lysis Buffer for SDS precipitation due to low storage temperature in which case it is necessary to dissolve the SDS by warming at 37 °C.

³ Before use, add 15 ml of absolute EtOH to the Washing buffer before use.

⁴ RNase A (10 mg/ml) can be stored for one year at room temperature (15-25 °C).

STORAGE CONDITION

The Fast DNA-spin™ Plasmid DNA Purification Kit can be stored at room temperature (15-25 °C) for up to 12 months without showing any reduction in performance and quality. For longer storage, this kit can be stored at 2-8 °C. If any precipitate forms in the buffers after storage at 2 - 8 °C, it should be dissolved by warming the buffers at 37 °C before use. After addition of RNase A and Lysis Viewer, Resuspension Buffer is stable for 6 months at 2 - 8 °C.

IMPORTANT NOTES

1. Add the provided RNase A and Lysis Viewer solutions to Resuspension Buffer before use, mix and store at 2 - 8 °C.
2. Check Lysis Buffer and Neutralization Buffer before use for salt precipitation. Redissolve any precipitate by warming at 37 °C.
3. Avoid direct contact of Lysis Buffer and Neutralization Buffer, immediately close the lid after use.
4. All centrifugation steps are carried out at 12,000rpm (~13,400 x g) in table-top microcentrifuge at room temperature (15 - 25 °C).
5. The obtained plasmid amount is influenced by bacteria culture density and plasmid copy number as well.
6. Lysis Viewer user guide : Lysis Viewer is an indicator which is harmless and used to make sure that the whole experimental process works well. Lysis Viewer should be

mixed with Resuspension Buffer in the ratio of 1 : 200 and the color of the mixed solution should be clear red. Add the mixed solution to cell culture and the solution would turn turbid red. After that, add Lysis Buffer to the turbid solution, the solution would turn clear purple which means a complete lysis. Add Neutralization Buffer to the purple solution and it would turn clear yellow, which indicate that the neutralization reaction has been done.

BACTERIAL CULTURE and COLLECTION

1. Inoculate 1 - 10ml of LB medium containing the appropriate antibiotic either with a single colony of transformed bacteria or with 0.2ml ~ 1.0ml of a small-scale LB culture grown from a single colony.
2. Incubate the culture at 37 °C with shaking for 12 - 16 hours.
 - If you use ampicillin as an antibiotic for culture (OD₆₀₀ 1.5 ~ 2.0), we recommend to increase your working ampicillin concentration up to 200 ~ 300 µg/mL to sustain selective antibiotic pressure for obtaining higher plasmid yield.
 - Growth for more than 16h is not recommended since cells begin to lyse and plasmid yields may be reduced. Use a tube or flask with a volume of at least 4 times volume of the culture.

PROTOCOL

1. Harvest 1 - 4ml bacterial cells in a microcentrifuge tube at 12,000rpm (~13,400 x g) in a conventional, table-top microcentrifuge for 1 min at room temperature (15 - 25 °C), then remove all traces of supernatant by inverting the open centrifuge tube until all medium has been drained.
2. Resuspend pelleted bacterial cells in 150 µl Resuspension Buffer by pipetting or vortex (Ensure that RNase A and Lysis Viewer have been added to Resuspension Buffer).

Note : Cell clumps indicate incomplete lysis, will result in lower yield and purity. Addition of Lysis Viewer will not have negative impact on following PCR, enzyme digestion and sequencing.
3. Add 150 µl Lysis Buffer and mix gently by inverting the tube 6 - 8 times.

Note : Mix gently by inverting the tube. Do not vortex, as this will result in shearing of genomic DNA. If necessary, continue inverting the tube until the solution becomes viscous and slightly clear. If not clear, probably due to incomplete lysis, please reduce the cells.
4. Add 350 µl Neutralization Buffer and mix immediately and quickly by inverting 12 - 20 times. The solution should become cloudy. Centrifuge for 2 min at 12,000rpm (~13,400 x g) in a table-top microcentrifuge.

Note : To avoid localized precipitation, mix the solution quickly, immediately after addition of Neutralization Buffer. The solution should be centrifuged again if there is still a lot of white precipitate can be seen in the supernatant. Since Lysis Viewer is applied, after the addition and mix of Neutralization Buffer, the solution should turn clear yellow. If there is still some purple liquid can be seen in the tube, keep inverting the tube until the color of solution turns completely clear yellow.
5. Transfer the supernatant from step 4 to the Spin Column (put in a Collection Tube) by pipetting. Centrifuge for 30 s at 12,000rpm (~13,400 x g). Discard the flow-through and set the Spin Column back into the Collection Tube.
6. Wash the Spin Column by adding 300 µl Washing Buffer (ensure the absolute ethanol has been added to Washing Buffer) and centrifuging for 30s at 12,000 rpm (~13,400 x g). Discard the flow-through, and put the Spin Column back into the Collection Tube.
7. Centrifuge for an additional 1 min at 12,000 rpm (~13,400 x g) to remove residual Washing Buffer.
8. Place the Spin Column in a clean 1.5 ml microcentrifuge tube. To elute DNA, add 50 - 100 µl Elution Buffer to the center of the Spin Column, centrifuge for 30 s at 12,000rpm (~13,400 x g).

Note : the volume of eluted buffer should not be less than 50 µl, otherwise it may affect recovery efficiency. The pH value of eluted buffer will have a great effect on eluting.

