

SiZer™ DNA Markers

DESCRIPTION

iNtRON supplies a wide range of products for accurate size and mass estimations (quantitation) of nucleic acid fragments. Nucleic Acid Markers are available for sizing linear, or supercoiled DNA and single-stranded RNA fragments. A variety of these markers are available in the ready-to-use SiZer™ formats.

SiZer™ DNA Markers are ideal for determining the size of double-stranded DNA from 60~10,000bp base pairs. The SiZer™ DNA Markers consist of 7 ~ 15linear double-stranded DNA fragments. Several fragments are present at increased intensity to allow easy identification. All fragments are precisely quantified and mixed during the production.

For 5 µl loading, all fragments except typical band DNA fragments are 40 ng. The typical band of DNA fragments is 100 ng. These ladders are pre-mixed with loading dye and are ready to use.

All DNA Markers can be stained with RedSafe™ Nucleic Acid Staining Solution, ethidium bromide (EtBr) or other DNA stains.

CHARACTERISTICS

- Ideal for determining the size of DNA
- Stable for more than 12 months at - 20 °C
- Ready to use without any handlings.

KIT CONTENTS

Product	Contents	Cat. No.
SiZer™-50 plus DNA Marker	0.5 ml	24072
SiZer™-100 DNA Marker	0.5 ml	24073
SiZer™-1000 DNA Marker	0.5 ml	24074
SiZer™-1000 plus DNA Marker	0.5 ml	24075

STORAGE

- Store at 4 °C and stable for more than 6 months. For more stable use, should be aliquoted and then stored at -20 °C. (stable for more than 12 months)

GENERAL USE

- No DNase and RNase detected.
- Load 5 µl per each well of Agarose gel.

QUALITY CONTROL

Well-defined bands are formed during agarose gel electrophoresis. The DNA concentration is determined spectrophotometrically. The absence of nucleases is confirmed by a direct nuclease activity assay.

ELECTROPHORESIS

- The 5 µl of ladder DNA was loaded, and then electrophoresed for 1hr at appropriate concentration of gel

PRODUCT USE LIMITATION

This product is developed, designed and sold exclusively for research purposes and in vitro use only. The product was not tested for use in diagnostics or for drug development, nor is it suitable for administration to humans or animals.

NOTICE BEFORE USE

- Do not heat before loading
- For quantification, adjust the concentration of the sample to equalize, it approximately with the amount of DNA in the nearest band of the ladder.
- Visualize DNA by staining RedSafe™, ethidium bromide (EtBr) or other DNA stains.


DETAIL INFORMATION

	Size range (bp)	Conc. (ng/µl)	Typical bands	Other bands	Loading Vol.	Band number	Contents
SiZer™ -50plus	50-500	128	100ng/5µl	40ng/5µl	5µl	13	50,100,150,200,250,300,400,500,600,700,800,900,1000
SiZer™ -100	100-1500	100	100ng/5µl	40ng/5µl	5µl	11	100,200,300,400,500,600,700,800,900,1000,1500
SiZer™ -1000	250-10000	120	100ng/5µl	40ng/5µl	5µl	12	250,500,750,1000,1500,2000,3000,4000,5000,6000,8000,10000
SiZer™ -1000 plus	100-10000	144	100ng/5µl	40ng/5µl	5µl	15	100,200,300,400,500,700,1000,1500,2000,3000,4000,5000,6000,8000,10000


RELATED PRODUCTS

Product Name	Cat.No.
RedSafe™ Nucleic Acid Staining Solution (20,000x)	21411
DNA-spin™ Plasmid DNA Extraction Kit	17096/17097/17098
MEGAquick-spin™ Total Fragment DNA Purification Kit	17286 / 17287/17288
Maxime™ PCR PreMix (i-StarTaq)	25165
Maxime™ PCR PreMix (i-pfu)	25185


SiZer™-50 plus



SiZer™-100



SiZer™-1000



SiZer™-1000 plus

