

News & Notice


January 2024

iNtRON Biotechnology

Let's make a time to meet "Special" PCR material products !



iNtRON One-Pick Products



RealMOD™ Probe R² 2X qRT-PCR mix (with UDG)

RealMOD™ Probe R² 2X qRT-PCR mix (with UDG) is a ready-to-use product that contains all the necessary reagents for qRT-PCR with UDG(Uracil DNA Glycosylase).

Characteristic

- Less cross contamination
- Detection and Qualification of RNA target
- Gene-expression analysis
- Pathogene detection & Qualification

Technical Data

RealMOD™ Probe R² 2X qRT-PCR mix (with UDG)




Figure 1. Amplification of various RNA using RealMOD™ Probe R² 2X qRT-PCR mix (w/UDG)

- A : TGE virus
- B : PED virus
- C : RdRP gene of SARS-CoV-2 virus
- D : E gene of SARS-CoV-2 virus
- E : N gene of SARS-CoV-2 virus

it is iNtRON.



• If you have any questions, please send an email to us. tech@intrnbio.com